

产品名称: AMPA Receptor 2 Rabbit Polyclonal Antibody
产品货号: APRab00369



产品概述 (Summary)

产品名称 (Production Name)	AMPA Receptor 2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRIA2
别名 (Alternative Names)	GRIA2; GLUR2; Glutamate receptor 2; GluR-2; AMPA-selective glutamate receptor 2; GluR-B; GluR-K2; Glutamate receptor ionotropic; AMPA 2; GluA2
基因 ID (Gene ID)	2891
蛋白 ID (SwissProt ID)	P42262.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 99 kDa; Observed MW: 99 kDa

研究背景 (Background)

产品名称: AMPA Receptor 2 Rabbit Polyclonal Antibody
产品货号: APRab00369




Ionotropic glutamate receptor. L-glutamate acts as an excitatory neurotransmitter at many synapses in the central nervous system. Binding of the excitatory neurotransmitter L-glutamate induces a conformation change, leading to the opening of the cation channel, and thereby converts the chemical signal to an electrical impulse. The receptor then desensitizes rapidly and enters a transient inactive state, characterized by the presence of bound agonist.


研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of AMPA Receptor 2/GluA2 in rat brain lysates using AMPA Receptor 2 antibody.



Western blot analysis of Glutamate receptor 2 in mouse brain lysates using Glutamate receptor 2 antibody.

注意事项 (Note)

For research use only .