

产品名称: Filamin A Rabbit Polyclonal Antibody
产品货号: APRab00370



产品概述 (Summary)

产品名称 (Production Name)	Filamin A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FLNA
别名 (Alternative Names)	ABP-280; ABPX; actin binding protein 280; Alpha-filamin; filamin 1; filamin A; alpha; FLN1; FLNA; FMD; MNS; NHBP; Non-muscle filamin; OPD; OPD1
基因 ID (Gene ID)	2316
蛋白 ID (SwissProt ID)	P21333.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 281 kDa; Observed MW: 281 kDa

研究背景 (Background)

产品名称: Filamin A Rabbit Polyclonal Antibody
产品货号: AP Rab00370




Promotes orthogonal branching of actin filaments and links actin filaments to membrane glycoproteins. Anchors various transmembrane proteins to the actin cytoskeleton and serves as a scaffold for a wide range of cytoplasmic signaling proteins. Interaction with FLNA may allow neuroblast migration from the ventricular zone into the cortical plate.


研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Filamin A in rat brain lysates using Filamin A antibody.



Western blot analysis of Filamin A in HepG2, 3T3, Hela lysates using Filamin A antibody.

注意事项 (Note)

For research use only .