

产品名称: CD80 Rabbit Polyclonal Antibody
产品货号: APRab00372



产品概述 (Summary)

产品名称 (Production Name)	CD80 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% protective protein
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD80 CD80; CD28LG; CD28LG1; LAB7; T-lymphocyte activation antigen CD80;
别名 (Alternative Names)	Activation B7-1 antigen; BB1; CTLA-4 counter-receptor B7.1; B7; CD antigen CD80
基因 ID (Gene ID)	941
蛋白 ID (SwissProt ID)	P33681.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 33 kDa; Observed MW: 65 kDa

产品名称: CD80 Rabbit Polyclonal Antibody
产品货号: AP Rab00372




研究背景 (Background)

Expressed on activated B-cells, macrophages and dendritic cells. Involved in the costimulatory signal essential for T-lymphocyte activation. T-cell proliferation and cytokine production is induced by the binding of CD28 or CTLA-4 to this receptor.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD80 in Raji lysates using CD80 antibody.

注意事项 (Note)

For research use only .