

产品名称: Fas Ligand Rabbit Polyclonal Antibody
产品货号: APRab00379



产品概述 (Summary)

产品名称 (Production Name)	Fas Ligand Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FASLG FASLG; APT1LG1; CD95L; FASL; TNFSF6; Tumor necrosis factor ligand
别名 (Alternative Names)	superfamily member 6; Apoptosis antigen ligand; APTL; CD95 ligand; CD95-L; Fas antigen ligand; Fas ligand; FasL; CD antigen CD178
基因 ID (Gene ID)	356
蛋白 ID (SwissProt ID)	P48023.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 31 kDa; Observed MW: 31 kDa

产品名称: Fas Ligand Rabbit Polyclonal Antibody
产品货号: APRab00379




研究背景 (Background)

Cytokine that binds to TNFRSF6/FAS, a receptor that transduces the apoptotic signal into cells.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of Fas Ligand in various lysates using Fas Ligand antibody.

注意事项 (Note)

For research use only .