

产品名称: ZAP70 Rabbit Polyclonal Antibody
产品货号: AP Rab00390



产品概述 (Summary)

产品名称 (Production Name)	ZAP70 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZAP70
别名 (Alternative Names)	ZAP70; SRK; Tyrosine-protein kinase ZAP-70; 70 kDa zeta-chain associated protein; Syk-related tyrosine kinase
基因 ID (Gene ID)	7535
蛋白 ID (SwissProt ID)	P43403.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa; Observed MW: 70 kDa

研究背景 (Background)

产品名称: ZAP70 Rabbit Polyclonal Antibody
产品货号: AP Rab00390




The Syk family protein tyrosine kinase Zap-70 is expressed in T and NK cells and plays a critical role in mediating T cell activation in response to T cell receptor (TCR) engagement. Following TCR engagement, Zap-70 is rapidly phosphorylated on several tyrosine residues through autophosphorylation and transphosphorylation by the Src family tyrosine kinase Lck.

研究领域 (Research Area)

Immunology

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue, using ZAP70 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.

注意事项 (Note)

For research use only .