

产品名称: kappa Actin Rabbit Polyclonal Antibody
产品货号: APRab00393



产品概述 (Summary)

产品名称 (Production Name)	kappa Actin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	ACTBL2
别名 (Alternative Names)	ACT
基因 ID (Gene ID)	345651
蛋白 ID (SwissProt ID)	Q562R1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 42 kDa; Observed MW: 42 kDa

研究背景 (Background)

Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all

产品名称: kappa Actin Rabbit Polyclonal Antibody
产品货号: APRab00393




eukaryotic cells.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of kappa Actin in various lysates using kappa Actin antibody.



Western blot analysis of kappa Actin in the HUVEC lysates using ACTBL2 antibody.

注意事项 (Note)

For research use only .