

**产品名称:** Factor VIII Rabbit Polyclonal Antibody

**产品货号:** APRab00399

## 产品概述 (Summary)

|                   |                                                                                                |
|-------------------|------------------------------------------------------------------------------------------------|
| 描述(Description)   | Rabbit polyclonal Antibody                                                                     |
| 宿主(Host)          | Rabbit                                                                                         |
| 应用(Application)   | WB,IHC,ELISA                                                                                   |
| 种属反应(Reactivity)  | Human,Mouse                                                                                    |
| 偶联物(Conjugation)  | Unconjugated                                                                                   |
| 修饰(Modification)  | Unmodified                                                                                     |
| 同种型(Isotype)      | IgG                                                                                            |
| 克隆(Clonality)     | Polyclonal                                                                                     |
| 剂型(Form)          | Liquid                                                                                         |
| 保存条件(Storage)     | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液(Buffer)      | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化 (Purification) | Affinity Purification                                                                          |

## 产品应用(Application)

|                       |                                                     |
|-----------------------|-----------------------------------------------------|
| 稀释比(Dilution Ratio)   | WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000 |
| 分子量(Molecular Weight) | Calculated MW: 267 kDa; Observed MW: 300 kDa        |

## 抗原信息(Antigen Information)

|                       |                                                                                                              |
|-----------------------|--------------------------------------------------------------------------------------------------------------|
| 基因名(Gene Name)        | F8                                                                                                           |
| 别名(Alternative Names) | F8; F8C; Coagulation factor VIII; Antihemophilic factor; AHF; Procoagulant component                         |
| 基因 ID(Gene ID)        | 2157                                                                                                         |
| SwissProt ID          | P00451                                                                                                       |
| 免疫原(Immunogen)        | The antiserum was produced against synthesized peptide derived from human Factor VIII.<br>AA range:2161-2210 |

## 研究背景 (Background)

F8: coagulation factor VIII, procoagulant component. This gene encodes coagulation factor VIII, which participates in the intrinsic pathway of blood coagulation; factor VIII is a cofactor for factor IXa which, in the presence of Ca<sup>2+</sup> and phospholipids,

converts factor X to the activated form Xa. This gene produces two alternatively spliced transcripts. Transcript variant 1 encodes a large glycoprotein, isoform a, which circulates in plasma and associates with von Willebrand factor in a noncovalent complex. This protein undergoes multiple cleavage events. Transcript variant 2 encodes a putative small protein, isoform b, which consists primarily of the phospholipid binding domain of factor VIIIc. This binding domain is essential for coagulant activity. Defects in this gene results in hemophilia A, a common recessive X-linked coagulation disorder.


### 研究领域 (Research Area)

Cardiovascular

### 注意事项 (Note)

For research use only.

### 图片 (Image Data)



Western blot analysis of Factor VIII in HuvEc lysates using Factor VIII antibody.

## EnkiLife 优势产品推荐

| 产品类别                        | 产品货号     | 产品名称                                               |
|-----------------------------|----------|----------------------------------------------------|
| WB 解决方案                     | RA10020  | 2h 极速 WB 即用型全流程试剂盒                                 |
|                             | RA10021  | 4h 快速 WB 即用型全流程试剂盒                                 |
|                             | RA10042  | 5h 畅享版 WB 全流程试剂盒                                   |
|                             | RA10037  | 校准级彩色预染蛋白 Marker (8-180kDa)                        |
|                             | RA10038  | 校准级彩色预染蛋白 Marker (10-250kDa)                       |
|                             | RA10039  | 校准级高分子彩色预染蛋白 Marker (25-400kDa)                    |
| TSA 多重荧光染色试剂盒               | RA10008  | TSA 双标三色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10009  | TSA 三标四色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10010  | TSA 四标五色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10011  | TSA 五标六色多重荧光染色试剂盒 (mIHC)                           |
|                             | RA10012  | TSA 六标七色多重荧光染色试剂盒 (mIHC)                           |
| IHC 检测试剂盒                   | RA10006  | HRP Anti-Mouse/Rabbit IHC Detection System         |
|                             | RA10007  | Polymer-HRP Anti-Mouse/Rabbit IHC Detection System |
| 抗体标记试剂盒                     | RE80004p | 辣根过氧化物酶(HRP)抗体标记试剂盒                                |
|                             | RE80002q | Sulfo-NHS-生物素标记试剂盒                                 |
|                             | RE80007p | Cy3 荧光素标记试剂盒                                       |
|                             | RE80011p | Fluor488 荧光素标记试剂盒                                  |
|                             | RE80017p | Fluor750 荧光素标记试剂盒                                  |
|                             | RE80005p | 藻红蛋白(R-PE) 抗体快速标记试剂盒                               |
|                             | RE80040  | PE-Cy7 串联染料抗体快速标记试剂盒                               |
| 稳转细胞系构建服务<br>(免费赠送全膜 WB 验证) | TS-0001  | 过表达稳转细胞系构建                                         |
|                             | TS-0002  | 敲低稳转细胞系构建                                          |
|                             | TS-0003  | 敲除细胞系构建                                            |

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>