

产品名称: Factor VIII Rabbit Polyclonal Antibody
产品货号: APRab00399



产品概述 (Summary)

产品名称 (Production Name)	Factor VIII Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	F8
别名 (Alternative Names)	F8; F8C; Coagulation factor VIII; Antihemophilic factor; AHF; Procoagulant component
基因 ID (Gene ID)	2157
蛋白 ID (SwissProt ID)	P00451.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 267 kDa; Observed MW: 300 kDa

研究背景 (Background)

产品名称: Factor VIII Rabbit Polyclonal Antibody
产品货号: APRab00399




F8: coagulation factor VIII, procoagulant component. This gene encodes coagulation factor VIII, which participates in the intrinsic pathway of blood coagulation; factor VIII is a cofactor for factor IXa which, in the presence of Ca^{+2} and phospholipids, converts factor X to the activated form Xa. This gene produces two alternatively spliced transcripts. Transcript variant 1 encodes a large glycoprotein, isoform a, which circulates in plasma and associates with von Willebrand factor in a noncovalent complex. This protein undergoes multiple cleavage events. Transcript variant 2 encodes a putative small protein, isoform b, which consists primarily of the phospholipid binding domain of factor VIIIc. This binding domain is essential for coagulant activity. Defects in this gene results in hemophilia A, a common recessive X-linked coagulation disorder.


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of Factor VIII in HuvEc lysates using Factor VIII antibody.



Western blot analysis of Factor VIII in HuvEc lysates using Factor VIII antibody.

注意事项 (Note)

For research use only .