

产品名称: GPR158 Rabbit Polyclonal Antibody
产品货号: APRab00401



产品概述 (Summary)

产品名称 (Production Name)	GPR158 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR158
别名 (Alternative Names)	Probable G-protein coupled receptor 158
基因 ID (Gene ID)	57512
蛋白 ID (SwissProt ID)	Q5T848.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 135 kDa; Observed MW: 150 kDa

研究背景 (Background)

Orphan receptor.


产品名称: GPR158 Rabbit Polyclonal Antibody
产品货号: AP Rab00401




研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Western blot analysis of GPR158 in COLO, 293T, HELA lysates using GPR158 antibody.




Immunofluorescence analysis of GPR158 in HUVEC using GPR158 antibody. The picture on the right is blocked using the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded Human brain tissue using GPR158 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.

产品名称: GPR158 Rabbit Polyclonal Antibody
产品货号: APRab00401



Western blot analysis of GPR158 in K562 lysates using GPR158 antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of GPR158 in HeLa using GPR158 antibody.

注意事项 (Note)

For research use only .