

产品名称: MKNK2 Rabbit Polyclonal Antibody
产品货号: APRab00407



产品概述 (Summary)

产品名称 (Production Name)	MKNK2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MKNK2 MKNK2; GPRK7; MNK2; MAP kinase-interacting serine/threonine-protein kinase 2; MAP kinase signal-integrating kinase 2; MAPK signal-integrating kinase 2; Mnk2
别名 (Alternative Names)	
基因 ID (Gene ID)	2872
蛋白 ID (SwissProt ID)	Q9HBH9.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 52 kDa

产品名称: MKNK2 Rabbit Polyclonal Antibody
产品货号: APRab00407




研究背景 (Background)

Serine/threonine-protein kinase that phosphorylates SFPQ/PSF, HNRNPA1 and EIF4E.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of MKNK2 in HeLa lysates using MKNK2 antibody.

注意事项 (Note)

For research use only .