

产品名称: Myomesin 2 Rabbit Polyclonal Antibody
产品货号: APRab00409



产品概述 (Summary)

产品名称 (Production Name)	Myomesin 2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYOM2
别名 (Alternative Names)	TTNAP
基因 ID (Gene ID)	9172
蛋白 ID (SwissProt ID)	P54296.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 165 kDa; Observed MW: 165 kDa

研究背景 (Background)

Major component of the vertebrate myofibrillar M band. Binds myosin, titin, and light meromyosin. This binding is dose

产品名称: Myomesin 2 Rabbit Polyclonal Antibody
产品货号: APRab00409




dependent.


研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Myomesin 2 in various lysates using Myomesin 2 antibody.



Western blot analysis of Myomesin 2 in HT-29 lysates using Myomesin2 antibody.



Western blot analysis of Myomesin 2 in COS7 and HT-29 lysates using MYOM2 antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Myomesin 2 Rabbit Polyclonal Antibody
产品货号: APRab00409



注意事项 (Note)

For research use only .