

产品名称: PKNOX2 Rabbit Polyclonal Antibody
产品货号: APRab00411



产品概述 (Summary)

产品名称 (Production Name)	PKNOX2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PKNOX2
别名 (Alternative Names)	PKNOX2; PREP2; Homeobox protein PKNOX2; Homeobox protein PREP-2; PBX/knotted homeobox 2
基因 ID (Gene ID)	63876
蛋白 ID (SwissProt ID)	Q96KN3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 52 kDa; Observed MW: 55 kDa

研究背景 (Background)


产品名称: PKNOX2 Rabbit Polyclonal Antibody
产品货号: APRab00411



Homeodomain proteins are sequence-specific transcription factors that share a highly conserved DNA-binding domain and play fundamental roles in cell proliferation, differentiation, and death.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of PKNOX2 in HepG2 lysates using PKNOX2 antibody.

注意事项 (Note)

For research use only .