

产品名称: Cytochrome P450 4Z1 Rabbit Polyclonal Antibody
产品货号: APRab00422



产品概述 (Summary)

产品名称 (Production Name)	Cytochrome P450 4Z1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYP4Z1
别名 (Alternative Names)	CYP4A20
基因 ID (Gene ID)	199974
蛋白 ID (SwissProt ID)	Q86W10.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

A cytochrome P450 monooxygenase that catalyzes the in-chain oxidation of fatty acids.


产品名称: Cytochrome P450 4Z1 Rabbit Polyclonal Antibody
产品货号: APRab00422




研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Immunofluorescence analysis of Cytochrome P450 4Z1 in HUVEC cells using Cytochrome P450 4Z1 antibody. Sample with blocking peptide on the right.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Cytochrome P450 4Z1 antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .