

产品名称: NCAPG2 Rabbit Polyclonal Antibody
产品货号: APRab00426



产品概述 (Summary)

产品名称 (Production Name)	NCAPG2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	NCAPG2
别名 (Alternative Names)	MTB; CAPG2; LUZP5; CAP-G2; hCAP-G2
基因 ID (Gene ID)	54892
蛋白 ID (SwissProt ID)	Q86XI2.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Regulatory subunit of the condensin-2 complex, a complex which establishes mitotic chromosome architecture and is

产品名称: NCAPG2 Rabbit Polyclonal Antibody
产品货号: APRab00426




involved in physical rigidity of the chromatid axis.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human lung carcinoma, using NCAPG2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.

注意事项 (Note)

For research use only .