

产品名称: Ubinuclein 1 Rabbit Polyclonal Antibody
产品货号: APRab00429



产品概述 (Summary)

产品名称 (Production Name)	Ubinuclein 1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Chromatography

免疫原信息 (Immunogen)

基因名 (Gene Name)	UBN1
别名 (Alternative Names)	VT; VT4
基因 ID (Gene ID)	29855
蛋白 ID (SwissProt ID)	Q9NPG3.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Acts as a novel regulator of senescence. Involved in the formation of senescence-associated heterochromatin foci (SAHF),

产品名称: Ubinuclein 1 Rabbit Polyclonal Antibody
产品货号: APRab00429




which represses expression of proliferation-promoting genes. Binds to proliferation-promoting genes. May be required for replication-independent chromatin assembly.


研究领域 (Research Area)

Microbiology

图片 (Image Data)



Immunofluorescence analysis of Ubinuclein 1 in A549 cells using Ubinuclein 1 antibody. Sample with blocking peptide on the right.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Ubinuclein 1 antibody. Sample with blocking peptide on the right. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .