

产品名称: GAS41 Rabbit Polyclonal Antibody
产品货号: APRab00430



产品概述 (Summary)

产品名称 (Production Name)	GAS41 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	YEATS4
别名 (Alternative Names)	YEATS4; GAS41; YEATS domain-containing protein 4; Glioma-amplified sequence 41; Gas41; NuMA-binding protein 1; NuBI-1; NuBI1
基因 ID (Gene ID)	8089
蛋白 ID (SwissProt ID)	O95619.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 26 kDa; Observed MW: 26 kDa

研究背景 (Background)


产品名称: GAS41 Rabbit Polyclonal Antibody
产品货号: APRab00430



YEATS domain containing 4 (YEATS4) Homo sapiens The protein encoded by this gene is found in the nucleoli. It has high sequence homology to human MLLT1, and yeast and human MLLT3 proteins. Both MLLT1 and MLLT3 proteins belong to a class of transcription factors, indicating that the encoded protein might also represent a transcription factor. This protein is thought to be required for RNA transcription. This gene has been shown to be amplified in tumors. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2014],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of GAS41 in HT-29 lysates using GAS41 antibody.

注意事项 (Note)

For research use only .