

产品名称: NT5DC3 Rabbit Polyclonal Antibody
产品货号: APRab00433



产品概述 (Summary)

产品名称 (Production Name)	NT5DC3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NT5DC3 NT5C3; P5N1; UMPH1; HSPC233; Cytosolic 5'-nucleotidase 3; Cytosolic 5'-nucleotidase III; cN-III; Pyrimidine 5'-nucleotidase 1; P5'N-1; P5N-1; PN-I; Uridine 5'-monophosphate hydrolase 1; p36
别名 (Alternative Names)	
基因 ID (Gene ID)	51559
蛋白 ID (SwissProt ID)	Q86UY8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 63 kDa; Observed MW: 60 kDa

产品名称: NT5DC3 Rabbit Polyclonal Antibody
产品货号: AP Rab00433




研究背景 (Background)

Predicted to enable 5'-nucleotidase activity. Predicted to be involved in dephosphorylation. Part of receptor complex.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NT5DC3 in A431 lysates using NT5DC3 antibody.

注意事项 (Note)

For research use only .