

产品名称: Alliinase Rabbit Polyclonal Antibody
产品货号: APRab00434



产品概述 (Summary)

产品名称 (Production Name)	Alliinase Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Plants

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	-
别名 (Alternative Names)	Alliin lyase 1
基因 ID (Gene ID)	-
蛋白 ID (SwissProt ID)	Q01594.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 56 kDa; Observed MW: 50 kDa

研究背景 (Background)

Able to cleave the C-S bond of sulfoxide derivatives of Cys to produce allicin, thus giving rise to all sulfur compounds which

产品名称: Alliinase Rabbit Polyclonal Antibody
产品货号: APRab00434




are responsible for most of the properties of garlic, such as the specific smell and flavor as well as the health benefits like blood lipid or blood pressure lowering.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Alliinase in purified alliinase lysates using Alliinase antibody.

注意事项 (Note)

For research use only .