

产品名称: MG53 Rabbit Polyclonal Antibody
产品货号: APRab00440



产品概述 (Summary)

产品名称 (Production Name)	MG53 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRIM72
别名 (Alternative Names)	Mitsugumin-53; TRI72; Tripartite motif containing 72; LOC493829
基因 ID (Gene ID)	493829
蛋白 ID (SwissProt ID)	Q6ZMU5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 53 kDa; Observed MW: 53 kDa

研究背景 (Background)

Muscle-specific protein that plays a central role in cell membrane repair by nucleating the assembly of the repair machinery

产品名称: MG53 Rabbit Polyclonal Antibody
产品货号: APRab00440




at injury sites.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of MG53 in mouse skeletal muscle lysates using MG53 antibody.



Immunohistochemistry analysis of paraffin-embedded mouse Pancreas Tissue using MG53 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .