

产品名称: Elongation Factor 1A2 Binding Protein Rabbit Polyclonal Antibody
产品货号: APRab00448



产品概述 (Summary)

产品名称 (Production Name)	Elongation Factor 1A2 Binding Protein Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IGFN1
别名 (Alternative Names)	IGFN1 antibody; immunoglobulin-like and fibronectin type III domain containing 1 antibody; EEF1A2BP1 antibody
基因 ID (Gene ID)	91156
蛋白 ID (SwissProt ID)	Q86VF2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 138 kDa; Observed MW: 115 kDa

产品名称: Elongation Factor 1A2 Binding Protein Rabbit
Polyclonal Antibody
产品货号: APRab00448


研究背景 (Background)

Predicted to be involved in homophilic cell adhesion via plasma membrane adhesion molecules; retina layer formation; and synapse assembly. Predicted to be located in Z disc and nucleus. Predicted to be active in synapse.

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of Elongation Factor 1A2 Binding Protein in Pig Skeletal Muscle, rat Skeletal Muscle, rat Skeletal Muscle lysates using eEF1A2 antibody.

注意事项 (Note)

For research use only .