

产品名称: GAL4 DNA Binding Domain Rabbit Polyclonal Antibody
产品货号: APRab00453



产品概述 (Summary)

产品名称 (Production Name)	GAL4 DNA Binding Domain Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Saccharomyces Cerevisiae

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GAL4
别名 (Alternative Names)	YPL248C
基因 ID (Gene ID)	855828.0
蛋白 ID (SwissProt ID)	P04386.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Refer to figures

研究背景 (Background)

产品名称: GAL4 DNA Binding Domain Rabbit Polyclonal Antibody
产品货号: APRab00453



This protein is a positive regulator for the gene expression of the galactose-induced genes such as GAL1, GAL2, GAL7, GAL10, and MEL1 which code for the enzymes used to convert galactose to glucose. It recognizes a 17 base pair sequence in (5'-CGGRNNRCYNYNCNCCG-3') the upstream activating sequence (UAS-G) of these genes.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of GAL4 DNA Binding Domain in Recombinant BD protein lysates using GAL4 DNA Binding Domain antibody.

注意事项 (Note)

For research use only .