

产品名称: GPR75 Rabbit Polyclonal Antibody  
产品货号: APRab00456



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | GPR75 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB, ICC/IF, ELISA                |
| 种属反应性 (Reactivity)     | Human, Mouse                     |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                        |
|------------------------|----------------------------------------|
| 基因名 (Gene Name)        | GPR75                                  |
| 别名 (Alternative Names) | Probable G-protein coupled receptor 75 |
| 基因 ID (Gene ID)        | 10936                                  |
| 蛋白 ID (SwissProt ID)   | O95800.                                |

## 产品应用 (Application)

|                          |                                                          |
|--------------------------|----------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 59 kDa; Observed MW: 59 kDa               |

## 研究背景 (Background)

G protein-coupled receptor that is activated by the chemokine CCL5/RANTES.


产品名称: GPR75 Rabbit Polyclonal Antibody  
产品货号: APRab00456




## 研究领域 (Research Area)

Signal Transduction

## 图片 (Image Data)



Western blot analysis of GPR75 in HT-29 and K562 lysates using GPR75 antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of GPR75 in LOVO using GPR75 antibody. The picture on the right is blocked using the synthesized peptide.

## 注意事项 (Note)

For research use only .