

产品名称: LIM Kinase 1 Rabbit Polyclonal Antibody
产品货号: APRab00458



产品概述 (Summary)

产品名称 (Production Name)	LIM Kinase 1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LIMK1
别名 (Alternative Names)	LIMK1; LIMK; LIM domain kinase 1; LIMK-1
基因 ID (Gene ID)	3984
蛋白 ID (SwissProt ID)	P53667.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW: 73 kDa; Observed MW: 73 kDa

研究背景 (Background)

Protein kinase which regulates actin filament dynamics. Phosphorylates and inactivates the actin binding/depolymerizing

产品名称: LIM Kinase 1 Rabbit Polyclonal Antibody
产品货号: APRab00458




factor cofilin, thereby stabilizing the actin cytoskeleton. Stimulates axonal outgrowth and may be involved in brain development. Isoform 3 has a dominant negative effect on actin cytoskeletal changes.

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of LIM Kinase 1 in mouse brain lysates using LIM Kinase 1 antibody.

注意事项 (Note)

For research use only .