

产品名称: MDM2 Rabbit Polyclonal Antibody
产品货号: APRab00459



产品概述 (Summary)

产品名称 (Production Name)	MDM2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MDM2
别名 (Alternative Names)	Double minute 2 protein; Hdm2; Oncoprotein Mdm2
基因 ID (Gene ID)	4193
蛋白 ID (SwissProt ID)	Q00987.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 90 kDa

研究背景 (Background)

E3 ubiquitin-protein ligase that mediates ubiquitination of p53/TP53, leading to its degradation by the proteasome.


产品名称: MDM2 Rabbit Polyclonal Antibody
产品货号: APRab00459




研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Western blot analysis of MDM2 in SKOV-3 lysates using MDM2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast carcinoma tissue using MDM2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval. Sample with blocking peptide on the right.



Immunofluorescence analysis of MDM2 in A549 cells using MDM2 antibody (green).

产品名称: MDM2 Rabbit Polyclonal Antibody
产品货号: APRab00459



注意事项 (Note)

For research use only .