

产品名称: CD241 Rabbit Polyclonal Antibody  
产品货号: APRab00470



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | CD241 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,IHC,ELISA                     |
| 种属反应性 (Reactivity)     | Human                            |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                               |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | RHAG<br>RHAG; RH50; Ammonium transporter Rh type A; Erythrocyte membrane glycoprotein Rh50; Erythrocyte plasma membrane 50 kDa glycoprotein;                                                                  |
| 别名 (Alternative Names) | Rh50A; Rhesus blood group family type A glycoprotein; Rh family type A glycoprotein; Rh type A glycoprotein; Rhesus blood group-associated ammonia channel; Rhesus blood group-associated glycoprotein; CD241 |
| 基因 ID (Gene ID)        | 6005                                                                                                                                                                                                          |
| 蛋白 ID (SwissProt ID)   | Q02094.                                                                                                                                                                                                       |

## 产品应用 (Application)

|                      |                                                     |
|----------------------|-----------------------------------------------------|
| 稀释比 (Dilution Ratio) | WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000 |
|----------------------|-----------------------------------------------------|

产品名称: CD241 Rabbit Polyclonal Antibody  
产品货号: APRab00470



蛋白分子量 (Molecular Weight)      Calculated MW: 44 kDa; Observed MW: 44 kDa


## 研究背景 (Background)

Associated with rhesus blood group antigen expression.


## 研究领域 (Research Area)

Cardiovascular

## 图片 (Image Data)




Western blot analysis of CD241 in K562 lysates using CD241 antibody.




Immunohistochemistry analysis of paraffin-embedded Human liver using CD241 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.


产品名称: CD241 Rabbit Polyclonal Antibody  
产品货号: APRab00470



Immunohistochemistry analysis of paraffin-embedded Human liver using CD241 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human brain using CD241 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using CD241 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**产品名称: CD241 Rabbit Polyclonal Antibody**  
**产品货号: APRab00470**



Immunohistochemistry analysis of paraffin-embedded Human brain using CD241 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

#### **注意事项 (Note)**

For research use only .