

产品名称: Inhibin beta E Rabbit Polyclonal Antibody
产品货号: APRab00474



产品概述 (Summary)

产品名称 (Production Name)	Inhibin beta E Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	INHBE
别名 (Alternative Names)	INHBE; Inhibin beta E chain; Activin beta-E chain
基因 ID (Gene ID)	83729
蛋白 ID (SwissProt ID)	P58166.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 39 kDa; Observed MW: 39 kDa

研究背景 (Background)

Inhibins and activins inhibit and activate, respectively, the secretion of follitropin by the pituitary gland.


产品名称: Inhibin beta E Rabbit Polyclonal Antibody
产品货号: APRab00474




研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of Inhibin beta E in L929 lysates using Inhibin beta E antibody.




Immunohistochemistry analysis of paraffin-embedded Human liver using Inhibin beta E antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded Human liver using Inhibin beta E antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Inhibin beta E Rabbit Polyclonal Antibody
产品货号: APRab00474



Immunohistochemistry analysis of paraffin-embedded mouse muscle using Inhibin beta E antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse muscle using Inhibin beta E antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .