

产品名称: MYH Rabbit Polyclonal Antibody
产品货号: APRab00476



产品概述 (Summary)

产品名称 (Production Name)	MYH Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MYH1 MYH1; Myosin-1; Myosin heavy chain 1; Myosin heavy chain 2x; MyHC-2x; Myosin heavy chain IIx/d; MyHC-IIx/d; Myosin heavy chain; skeletal muscle; adult 1; MYH2; MYHSA2; Myosin-2; Myosin heavy chain 2Myosin heavy chain 2; Myosin heavy chain 2a; MyHC-2a; Myosin heavy chain IIa; MyHC-IIa; Myosin heavy chain; skeletal muscle; adult 2; MYH3; Myosin-3; Muscle embryonic myosin heavy chain; Myosin heavy chain 3
别名 (Alternative Names)	
基因 ID (Gene ID)	4619/4620/4621/4622/4624/4625/4626
蛋白 ID (SwissProt ID)	P12882/Q9UKX2/P11055/Q9Y623/P13533/P12883/P13535.

产品应用 (Application)

产品名称: MYH Rabbit Polyclonal Antibody
产品货号: APRab00476



稀释比 (Dilution Ratio)

WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000

蛋白分子量 (Molecular Weight)

Calculated MW: 223 kDa; Observed MW: 223 kDa


研究背景 (Background)

Muscle contraction.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Western blot analysis of MYH pAb in 4T1 lysates using MYH antibody.




Immunohistochemistry analysis of paraffin-embedded Human muscle using MYH antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.


产品名称: MYH Rabbit Polyclonal Antibody
产品货号: APRab00476



Immunohistochemistry analysis of paraffin-embedded Human muscle using MYH antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse muscle using MYH antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse muscle using MYH antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .