

产品名称: FGL1 Rabbit Polyclonal Antibody
产品货号: APRab00485



产品概述 (Summary)

产品名称 (Production Name)	FGL1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGL1 FGL1; HFREP1; Fibrinogen-like protein 1; HP-041; Hepassocin; Hepatocyte-derived fibrinogen-related protein 1; HFREP-1; Liver fibrinogen-related protein 1; LFIRE-1
别名 (Alternative Names)	
基因 ID (Gene ID)	2267
蛋白 ID (SwissProt ID)	Q08830.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 36 kDa; Observed MW: 36 kDa

产品名称: FGL1 Rabbit Polyclonal Antibody
产品货号: APRab00485



研究背景 (Background)

Immune suppressive molecule that inhibits antigen-specific T-cell activation by acting as a major ligand of LAG3.


研究领域 (Research Area)

Cardiovascular


图片 (Image Data)



Western blot analysis of FGL1 in mouse liver lysates using FGL1 antibody.



Immunohistochemistry analysis of paraffin-embedded rat liver using FGL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemistry analysis of paraffin-embedded rat liver using FGL1 antibody. High-pressure and temperature Sodium

产品名称: FGL1 Rabbit Polyclonal Antibody
产品货号: APRab00485



Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded mouse liver using FGL1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .