

产品名称: Apelin Rabbit Polyclonal Antibody
产品货号: APRab00487



产品概述 (Summary)

产品名称 (Production Name)	Apelin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APLN
别名 (Alternative Names)	APLN; APEL; Apelin; APJ endogenous ligand
基因 ID (Gene ID)	8862
蛋白 ID (SwissProt ID)	Q9ULZ1.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Endogenous ligand for the apelin receptor (APLNR).


产品名称: Apelin Rabbit Polyclonal Antibody
产品货号: AP Rab00487




研究领域 (Research Area)

Neuroscience


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded rat brain using Apelin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Apelin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded rat brain using Apelin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Apelin Rabbit Polyclonal Antibody
产品货号: APRab00487



注意事项 (Note)

For research use only .