

产品名称: PPP1R15B Rabbit Polyclonal Antibody
产品货号: APRab00501



产品概述 (Summary)

产品名称 (Production Name)	PPP1R15B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP1R15B
别名 (Alternative Names)	PPP1R15B; Protein phosphatase 1 regulatory subunit 15B
基因 ID (Gene ID)	84919.0
蛋白 ID (SwissProt ID)	Q5SWA1.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 79 kDa; Observed MW: 79 kDa

研究背景 (Background)

Maintains low levels of EIF2S1 phosphorylation in unstressed cells by promoting its dephosphorylation by PP1.


产品名称: PPP1R15B Rabbit Polyclonal Antibody
产品货号: APRab00501




研究领域 (Research Area)

Signal Transduction


图片 (Image Data)



Western blot analysis of PPP1R15B in NIH3T3 lysates using PPP1R15B antibody.




Immunohistochemistry analysis of paraffin-embedded Human liver using PPP1R15B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using PPP1R15B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: PPP1R15B Rabbit Polyclonal Antibody
产品货号: APRab00501



Immunohistochemistry analysis of paraffin-embedded Human stomach using PPP1R15B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .