

产品名称: MIA2 Rabbit Polyclonal Antibody
产品货号: APRab00506



产品概述 (Summary)

产品名称 (Production Name)	MIA2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MIA2
别名 (Alternative Names)	Melanoma inhibitory activity protein 2
基因 ID (Gene ID)	4253
蛋白 ID (SwissProt ID)	Q96PC5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 160 kDa; Observed MW: 65 kDa

研究背景 (Background)

Plays a role in the transport of cargos that are too large to fit into COPII-coated vesicles and require specific mechanisms to

产品名称: MIA2 Rabbit Polyclonal Antibody
产品货号: APRab00506




be incorporated into membrane-bound carriers and exported from the endoplasmic reticulum.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of MIA2 in HEPG2 lysates using MIA2 antibody.




Immunohistochemistry analysis of paraffin-embedded Human liver using MIA2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human pancreas using MIA2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: MIA2 Rabbit Polyclonal Antibody
产品货号: APRab00506



Immunohistochemistry analysis of paraffin-embedded Human pancreas using MIA2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .