

产品名称: Growth Hormone Rabbit Polyclonal Antibody
产品货号: AP Rab00507



产品概述 (Summary)

产品名称 (Production Name)	Growth Hormone Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GH1/2 Somatotropin (Growth hormone; GH; GH-N; Growth hormone 1; Pituitary growth hormone)
别名 (Alternative Names)	Growth hormone variant (GH-V; Growth hormone 2; Placenta-specific growth hormone)
基因 ID (Gene ID)	2689/2688
蛋白 ID (SwissProt ID)	P01241/P01242.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

产品名称: Growth Hormone Rabbit Polyclonal Antibody
产品货号: APRab00507




研究背景 (Background)

The protein encoded by this gene is a member of the somatotropin/prolactin family of hormones which play an important role in growth control.


研究领域 (Research Area)

Signal Transduction

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded Human placenta using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.




Immunohistochemical analysis of paraffin-embedded Human tonsils using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Growth Hormone Rabbit Polyclonal Antibody
产品货号: APRab00507



Immunohistochemistry analysis of paraffin-embedded Human brain using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human brain using Growth Hormone antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .