

产品名称: Integrin alpha 4 Rabbit Polyclonal Antibody
产品货号: APRab00523



产品概述 (Summary)

产品名称 (Production Name)	Integrin alpha 4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ITGAD
别名 (Alternative Names)	ITGAD; Integrin alpha-D; ADB2; CD11 antigen-like family member D; Leukointegrin alpha D; CD11d
基因 ID (Gene ID)	3681
蛋白 ID (SwissProt ID)	Q13349.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 127 kDa; Observed MW: 127 kDa

研究背景 (Background)

产品名称: Integrin alpha 4 Rabbit Polyclonal Antibody
产品货号: APRab00523




Integrin alpha-D/beta-2 is a receptor for ICAM3 and VCAM1.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Western blot analysis of Integrin alpha 4 in 4T 293 lysates using Integrin alpha 4 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung using Integrin alpha 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human spleen using Integrin alpha 4 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

产品名称: Integrin alpha 4 Rabbit Polyclonal Antibody
产品货号: APRab00523



注意事项 (Note)

For research use only .