

**产品名称: GDF6 Rabbit Polyclonal Antibody**  
**产品货号: APRab00530**



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | GDF6 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | WB,ELISA                        |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                 |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                                                       |
|------------------------|---------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | GDF6                                                                                  |
| 别名 (Alternative Names) | GDF6; GDF16; Growth/differentiation factor 6; GDF-6; Growth/differentiation factor 16 |
| 基因 ID (Gene ID)        | 392255                                                                                |
| 蛋白 ID (SwissProt ID)   | Q6KF10.                                                                               |

## 产品应用 (Application)

|                          |                                            |
|--------------------------|--------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,ELISA 1:5000-1:20000       |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 51 kDa; Observed MW: 51 kDa |

## 研究背景 (Background)

产品名称: GDF6 Rabbit Polyclonal Antibody  
产品货号: APRab00530




Growth factor that controls proliferation and cellular differentiation in the retina and bone formation.


## 研究领域 (Research Area)

Signal Transduction


## 图片 (Image Data)



Western blot analysis of GDF6 in rat kidney lysates using GDF6 antibody.



Western blot analysis of GDF6 in rat kidney lysates using GDF6 antibody.



Western blot analysis of GDF6 in 293T lysates using GDF6 antibody.

**产品名称: GDF6 Rabbit Polyclonal Antibody**  
**产品货号: APRab00530**



---

**注意事项 (Note)**

For research use only .