

**产品名称: Dectin-2 Rabbit Polyclonal Antibody**  
**产品货号: AP Rab00532**



## 产品概述 (Summary)

|                        |                                     |
|------------------------|-------------------------------------|
| 产品名称 (Production Name) | Dectin-2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody          |
| 宿主 (Host)              | Rabbit                              |
| 应用 (Application)       | WB,IHC,ELISA                        |
| 种属反应性 (Reactivity)     | Human                               |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                                |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CLEC6A<br>Dectin-2; CLEC6A; CLECSF10; DECTIN2; C-type lectin domain family 6 member A; C-type lectin superfamily member 10; Dendritic cell-associated C-type lectin 2; DC-associated C-type lectin 2; Dectin-2 |
| 别名 (Alternative Names) |                                                                                                                                                                                                                |
| 基因 ID (Gene ID)        | 93978                                                                                                                                                                                                          |
| 蛋白 ID (SwissProt ID)   | Q6EIG7.                                                                                                                                                                                                        |

## 产品应用 (Application)

|                          |                                                     |
|--------------------------|-----------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 24 kDa; Observed MW: 24 kDa          |

产品名称: Dectin-2 Rabbit Polyclonal Antibody  
产品货号: APRab00532




## 研究背景 (Background)

Like dectin-1 (CLEC7A; MIM 606264), CLEC6A, or dectin-2, is a type II membrane receptor with an extracellular C-type lectin-like domain fold. However, unlike dectin-1, dectin-2 lacks an immunoreceptor tyrosine-based activation motif (ITAM) in its cytoplasmic domain


## 研究领域 (Research Area)

Immunology


## 图片 (Image Data)



Western blot analysis of Dectin-2 in SH-SY5Y, 293T lysates using Dectin-2 antibody.



Western blot analysis of Dectin-2 in K562 lysates using Dectin-2 antibody.



**产品名称: Dectin-2 Rabbit Polyclonal Antibody**  
**产品货号: APRab00532**



---

Immunohistochemistry analysis of paraffin-embedded Human heart using Dectin-2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

**注意事项 (Note)**

For research use only .