

产品名称: TP53INP2 Rabbit Polyclonal Antibody
产品货号: APRab00533



产品概述 (Summary)

产品名称 (Production Name)	TP53INP2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TP53INP2 TP53INP2; C20orf110; DOR; PINH; Tumor protein p53-inducible nuclear protein 2; Diabetes and obesity-regulated gene; p53-inducible protein U; PIG-U
别名 (Alternative Names)	
基因 ID (Gene ID)	58476
蛋白 ID (SwissProt ID)	Q8IXH6.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 24 kDa; Observed MW: 24 kDa

产品名称: TP53INP2 Rabbit Polyclonal Antibody
产品货号: APRab00533




研究背景 (Background)

Dual regulator of transcription and autophagy.


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of TP53INP2 in SH-SY5Y, 293T, mouse kidney lysates using TP53INP2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using TP53INP2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using TP53INP2 antibody. High-pressure and

产品名称: TP53INP2 Rabbit Polyclonal Antibody
产品货号: APRab00533



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .