

产品名称: CD66a/b/c Rabbit Polyclonal Antibody
产品货号: APRab00539



产品概述 (Summary)

产品名称 (Production Name)	CD66a/b/c Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEACAM1/CEACAM6/CEACAM8
别名 (Alternative Names)	BGP; BGP1; NCA; CGM6
基因 ID (Gene ID)	634/4680/1088
蛋白 ID (SwissProt ID)	P13688/P40199/P31997.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:100, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	-

研究背景 (Background)

Plays a role as coinhibitory receptor in immune response, insulin action and functions also as an activator during

产品名称: CD66a/b/c Rabbit Polyclonal Antibody
产品货号: APRab00539




angiogenesis.


研究领域 (Research Area)

Immunology


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded Human kidney using CD66a/b/c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using CD66a/b/c antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemistry analysis of paraffin-embedded Human stomach using CD66a/b/c antibody. High-pressure and

产品名称: CD66a/b/c Rabbit Polyclonal Antibody
产品货号: APRab00539



temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .