

产品名称: B4GALT5 Rabbit Polyclonal Antibody
产品货号: APRab00548



产品概述 (Summary)

产品名称 (Production Name)	B4GALT5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	B4GALT5 B4GALT5; Beta-1; 4-galactosyltransferase 5; Beta-1; 4-GalTase 5; Beta4Gal-T5;
别名 (Alternative Names)	b4Gal-T5; Beta-1; 4-GalT II; UDP-Gal:beta-GlcNAc beta-1; 4-galactosyltransferase 5; UDP-galactose:beta-N-acetylglucosamine beta-1; 4-galactosyltransferase 5
基因 ID (Gene ID)	9334
蛋白 ID (SwissProt ID)	O43286.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 40 kDa

产品名称: B4GALT5 Rabbit Polyclonal Antibody
产品货号: APRab00548




研究背景 (Background)

This gene is one of seven beta-1,4-galactosyltransferase (beta4GalT) genes. They encode type II membrane-bound glycoproteins that appear to have exclusive specificity for the donor substrate UDP-galactose; all transfer galactose in a beta1,4 linkage to similar acceptor sugars: GlcNAc, Glc, and Xyl. Each beta4GalT has a distinct function in the biosynthesis of different glycoconjugates and saccharide structures. As type II membrane proteins, they have an N-terminal hydrophobic signal sequence that directs the protein to the Golgi apparatus and which then remains uncleaved to function as a transmembrane anchor. By sequence similarity, the beta4GalTs form four groups: beta4GalT1 and beta4GalT2, beta4GalT3 and beta4GalT4, beta4GalT5 and beta4GalT6, and beta4GalT7. The function of the enzyme encoded by this gene is not clear. This gene was previously designated as B4GALT4 but was renamed to B4GALT5. In the literature it is also referred to as beta4GalT2.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of B4GALT5 in various lysates using B4GALT5 antibody.



Western blot analysis of B4GALT5 in 293 lysates using β1, 4GalT5 antibody.

产品名称: B4GALT5 Rabbit Polyclonal Antibody
产品货号: APRab00548



注意事项 (Note)

For research use only .