

产品名称: NRG2 Rabbit Polyclonal Antibody
产品货号: APRab00553



产品概述 (Summary)

产品名称 (Production Name)	NRG2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NRG2
别名 (Alternative Names)	NRG2; NTAK; Pro-neuregulin-2; membrane-bound isoform; Pro-NRG2
基因 ID (Gene ID)	9542
蛋白 ID (SwissProt ID)	O14511.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 92 kDa; Observed MW: 92 kDa

研究背景 (Background)

Direct ligand for ERBB3 and ERBB4 tyrosine kinase receptors.


产品名称: NRG2 Rabbit Polyclonal Antibody
产品货号: APRab00553




研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western blot analysis of NRG2 in HT-29 lysates using NRG2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human brain using NRG2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .