

产品名称: Placental Alkaline Phosphatase Rabbit Polyclonal Antibody
产品货号: APRab00558



产品概述 (Summary)

产品名称 (Production Name)	Placental Alkaline Phosphatase Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALPP/ALPPL2 ALPP; PLAP; Alkaline phosphatase; placental type; Alkaline phosphatase Regan isozyme; Placental alkaline phosphatase 1; PLAP-1; ALPPL2; ALPPL; Alkaline phosphatase; placental-likeALP-1; Alkaline phosphatase Nagao isozyme; Germ cell alkaline phosphatase; GCAP; Placental alkaline phosphatase-like; PLAP-like
别名 (Alternative Names)	
基因 ID (Gene ID)	250/251
蛋白 ID (SwissProt ID)	P05187/P10696.

产品应用 (Application)

产品名称: Placental Alkaline Phosphatase Rabbit Polyclonal Antibody
产品货号: APRab00558



稀释比 (Dilution Ratio) WB 1:500-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight) Calculated MW: 57 kDa; Observed MW: 60 kDa


研究背景 (Background)

PLAP may assist in guiding migratory cells and transporting specific molecules, such as fatty acids and immunoglobulins, across the plasma membrane. The three tissue-specific APs identified in human, PLAP, germ cell AP (GCAP) and intestinal AP, are 90-98% homologous and their genes are clustered on chromosome 2q.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of Placental Alkaline Phosphatase in 293 lysates using Placental Alkaline Phosphatase antibody.



Western blot analysis of Placental Alkaline Phosphatase in lysates using ALPP/ALPPL2 antibody.

注意事项 (Note)

For research use only .