

产品名称: CD42d Rabbit Polyclonal Antibody
产品货号: APRab00563



产品概述 (Summary)

产品名称 (Production Name)	CD42d Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GP5
别名 (Alternative Names)	GP5; Platelet glycoprotein V; GPV; Glycoprotein 5; CD42d
基因 ID (Gene ID)	2814
蛋白 ID (SwissProt ID)	P40197.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,IHC 1:50-1:100,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 61 kDa; Observed MW: 61 kDa

研究背景 (Background)

The GPIb-V-IX complex functions as the vWF receptor and mediates vWF-dependent platelet adhesion to blood vessels.


产品名称: CD42d Rabbit Polyclonal Antibody
产品货号: AP Rab00563




研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Western blot analysis of CD42d in HeLa lysates using CD42d antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsils using CD42d antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

注意事项 (Note)

For research use only .