

产品名称: CEACAM6 Rabbit Polyclonal Antibody
产品货号: APRab00564



产品概述 (Summary)

产品名称 (Production Name)	CEACAM6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEACAM3/CEACAM6 CEACAM3; CD66D; CGM1; Carcinoembryonic antigen-related cell adhesion molecule 3; Carcinoembryonic antigen CGM1; CD66d; CEACAM6; NCA;
别名 (Alternative Names)	Carcinoembryonic antigen-related cell adhesion molecule 6; Non-specific crossreacting antigen; Normal cross-reacting antigen; CD66c
基因 ID (Gene ID)	1084/4680
蛋白 ID (SwissProt ID)	P40198/P40199.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 27,37 kDa; Observed MW: 30 kDa

产品名称: CEACAM6 Rabbit Polyclonal Antibody
产品货号: APRab00564




研究背景 (Background)

This gene encodes a member of the family of carcinoembryonic antigen-related cell adhesion molecules (CEACAMs), which are used by several bacterial pathogens to bind and invade host cells. The encoded transmembrane protein directs phagocytosis of several bacterial species that is dependent on the small GTPase Rac. It is thought to serve an important role in controlling human-specific pathogens by the innate immune system. Alternatively spliced transcript variants have been described.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CEACAM6 in HepG2, mouse brain lysates using CEACAM6 antibody.



Western blot analysis of CEACAM6 in HEPG2 lysates using CD66c/d antibody.

注意事项 (Note)

For research use only .