

**产品名称: CD232 Rabbit Polyclonal Antibody**  
**产品货号: APRab00568**



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | CD232 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,ELISA                         |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                  |

## 产品性能 (Performance)

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                   |
| 修饰 (Modification)   | Unmodified                                                                                     |
| 同种型 (Isotype)       | IgG                                                                                            |
| 克隆 (Clonality)      | Polyclonal                                                                                     |
| 形式 (Form)           | Liquid                                                                                         |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.       |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3. |
| 纯化方式 (Purification) | Affinity Purification                                                                          |

## 免疫原信息 (Immunogen)

|                        |                                                                               |
|------------------------|-------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | PLXNC1                                                                        |
| 别名 (Alternative Names) | PLXNC1; VESPR; Plexin-C1; Virus-encoded semaphorin protein receptor;<br>CD232 |
| 基因 ID (Gene ID)        | 10154                                                                         |
| 蛋白 ID (SwissProt ID)   | O60486.                                                                       |

## 产品应用 (Application)

|                          |                                              |
|--------------------------|----------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:1000,ELISA 1:5000-1:20000         |
| 蛋白分子量 (Molecular Weight) | Calculated MW: 176 kDa; Observed MW: 176 kDa |

## 研究背景 (Background)

产品名称: CD232 Rabbit Polyclonal Antibody  
产品货号: APRab00568




Receptor for SEMA7A, for smallpox semaphorin A39R, vaccinia virus semaphorin A39R and for herpesvirus Sema protein.


## 研究领域 (Research Area)

Cancer

## 图片 (Image Data)



Western blot analysis of CD232 in 3T3 lysates using CD232 antibody.



Western blot analysis of CD232 in 3T3 lysates using CD232 antibody.

## 注意事项 (Note)

For research use only .