

产品名称: CD155 Rabbit Polyclonal Antibody
产品货号: APRab00569



产品概述 (Summary)

产品名称 (Production Name)	CD155 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PVR
别名 (Alternative Names)	PVR; PVS; Poliovirus receptor; Nectin-like protein 5; NECL-5; CD155
基因 ID (Gene ID)	5817
蛋白 ID (SwissProt ID)	P15151.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 45 kDa; Observed MW: 45 kDa

研究背景 (Background)

Mediates NK cell adhesion and triggers NK cell effector functions. Binds two different NK cell receptors: CD96 and CD226.

产品名称: CD155 Rabbit Polyclonal Antibody
产品货号: APRab00569




These interactions accumulates at the cell-cell contact site, leading to the formation of a mature immunological synapse between NK cell and target cell.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD155 in KB lysates using CD155 antibody.



Western blot analysis of CD155 in KB lysates using CD155 antibody.

注意事项 (Note)

For research use only .