

产品名称: CD300d Rabbit Polyclonal Antibody
产品货号: APRab00570



产品概述 (Summary)

产品名称 (Production Name)	CD300d Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD300LD
别名 (Alternative Names)	CD300LD; CD300D; CMRF35A4; CMRF35-like molecule 4; CLM-4; CD300 antigen-like family member D; CMRF35-A4; CD300d
基因 ID (Gene ID)	100131439
蛋白 ID (SwissProt ID)	Q6UXZ3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 22 kDa

研究背景 (Background)

产品名称: CD300d Rabbit Polyclonal Antibody
产品货号: APRab00570




Predicted to enable transmembrane signaling receptor activity and virus receptor activity. Predicted to be involved in immune system process. Predicted to act upstream of or within regulation of interleukin-6 production and regulation of tumor necrosis factor production. Predicted to be integral component of membrane. Predicted to be active in plasma membrane.


研究领域 (Research Area)

Immunology

图片 (Image Data)



Western blot analysis of CD300d in NIH3T3 lysates using CD300d antibody.



Western blot analysis of CD300d in 3T3 lysates using CD300d antibody.

注意事项 (Note)

For research use only .