

产品名称: MRPL35 Rabbit Polyclonal Antibody
产品货号: APRab00574



产品概述 (Summary)

产品名称 (Production Name)	MRPL35 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRPL35
别名 (Alternative Names)	MRPL35; BM-007; 39S ribosomal protein L35; mitochondrial; L35mt; MRP-L35
基因 ID (Gene ID)	51318
蛋白 ID (SwissProt ID)	Q9NZE8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 22 kDa; Observed MW: 22 kDa

研究背景 (Background)

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the

产品名称: MRPL35 Rabbit Polyclonal Antibody
产品货号: APRab00574




mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit.


研究领域 (Research Area)

Tags & Cell Markers

图片 (Image Data)



Western blot analysis of MRPL35 in 293T , COLO lysates using MRPL35 antibody.



Western blot analysis of MRPL35 in 293 lysates using MRPL35 antibody.

注意事项 (Note)

For research use only .