

产品名称: ApoL6 Rabbit Polyclonal Antibody
产品货号: APRab00600



产品概述 (Summary)

产品名称 (Production Name)	ApoL6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOL6
别名 (Alternative Names)	APOL6; Apolipoprotein L6; Apolipoprotein L-VI; ApoL-VI
基因 ID (Gene ID)	80830
蛋白 ID (SwissProt ID)	Q9BWW8.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 38 kDa; Observed MW: 38 kDa

研究背景 (Background)

This gene is a member of the apolipoprotein L gene family. The encoded protein is found in the cytoplasm, where it may

产品名称: ApoL6 Rabbit Polyclonal Antibody
产品货号: APRab00600




affect the movement of lipids or allow the binding of lipids to organelles. Widely expressed; highly expressed in the uterus, fetal brain and spinal cord, also detected in heart, liver, lung, colon, spleen, thymus, prostate, placenta, adrenal gland, salivary and mammary gland.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



Western blot analysis of ApoL6 in K562 lysates using ApoL6 antibody.

注意事项 (Note)

For research use only .