

产品名称: DNA Polymerase mu Rabbit Polyclonal Antibody
产品货号: APRab00606



产品概述 (Summary)

产品名称 (Production Name)	DNA Polymerase mu Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POLM
别名 (Alternative Names)	POLM; polmu; DNA-directed DNA/RNA polymerase mu; Pol Mu; Terminal transferase
基因 ID (Gene ID)	27434
蛋白 ID (SwissProt ID)	Q9NP87.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 55 kDa; Observed MW: 55 kDa

研究背景 (Background)

产品名称: DNA Polymerase mu Rabbit Polyclonal Antibody
产品货号: APRab00606




Gap-filling polymerase involved in repair of DNA double-strand breaks by non-homologous end joining (NHEJ).

研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Western blot analysis of DNA Polymerase mu in K562 lysates using DNA Pol μ antibody.

注意事项 (Note)

For research use only .