

产品名称: Wnt16 Rabbit Polyclonal Antibody
产品货号: APRab00607



产品概述 (Summary)

产品名称 (Production Name)	Wnt16 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	WNT16
别名 (Alternative Names)	WNT16; Protein Wnt-16
基因 ID (Gene ID)	51384
蛋白 ID (SwissProt ID)	Q9UBV4.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 41 kDa; Observed MW: 41 kDa

研究背景 (Background)

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a

产品名称: Wnt16 Rabbit Polyclonal Antibody
产品货号: APRab00607




signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters .

研究领域 (Research Area)

Stem Cells

图片 (Image Data)



Western blot analysis of Wnt16 in rat kidney lysates using Wnt16 antibody.

注意事项 (Note)

For research use only .