

产品名称: Thioredoxin Domain Containing Protein 17 Rabbit Polyclonal Antibody
产品货号: APRab00610



产品概述 (Summary)

产品名称 (Production Name)	Thioredoxin Domain Containing Protein 17 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% sodium azide, pH 7.3.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TXNDC17
别名 (Alternative Names)	TXNDC17; TXNL5; Thioredoxin domain-containing protein 17; 14 kDa thioredoxin-related protein; TRP14; Protein 42-9-9; Thioredoxin-like protein 5
基因 ID (Gene ID)	84817
蛋白 ID (SwissProt ID)	Q9BRA2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW: 14 kDa; Observed MW: 14 kDa

产品名称: Thioredoxin Domain Containing Protein 17 Rabbit Polyclonal Antibody
产品货号: APRab00610


研究背景 (Background)

Disulfide reductase. May participate in various redox reactions through the reversible oxidation of its active center dithiol to a disulfide and catalyze dithiol-disulfide exchange reactions. Modulates TNF-alpha signaling and NF-kappa-B activation. Has peroxidase activity and may contribute to the elimination of cellular hydrogen peroxide.

研究领域 (Research Area)

Signal Transduction

图片 (Image Data)



Western blot analysis of Thioredoxin Domain Containing Protein 17 in HL-60, Raji, K562, SW480, A549, 22RV1 lysates using Thioredoxin Domain Containing Protein 17 antibody.

注意事项 (Note)

For research use only .